

# Alaris® Syringe module-compatible syringes and flow rate ranges

## Alaris products

Medications that will be delivered on the Alaris® Syringe module must be prepared in a syringe that will accommodate the desired flow rate. The following chart shows the flow rate ranges for each syringe size. In addition, the compatible syringe sizes with re-order numbers are displayed.

Syringe size	Brand and model numbers				
	BD	Monoject	Terumo	IVAC®	Flow rate range
1 mL	309628	n/a	n/a	n/a	0.01 to 30 mL/hr
3 mL	309657	SoftPack only: 1180300777	SS-03L	n/a	0.01 to 100 mL/hr
5 mL	309646	n/a	SS-05L (6 mL total)	n/a	0.1 to 150 mL/hr
6 mL	n/a	SoftPack: 1180600777 Rigid Pack: 8881516937	SS-05L (6 mL total)	n/a	0.1 to 150 mL/hr
10 mL	309604* 300912	n/a	SS-10L (12 mL total)	n/a	0.1 to 250 mL/hr
12 mL	n/a	SoftPack: 1181200777 Rigid Pack: 8881512878	SS-10L (12 mL total)	n/a	0.1 to 250 mL/hr
20 mL	309661	SoftPack: 1182000777 Rigid Pack: 8881520657	SS-20L SS-20L2	n/a	0.1 to 500 mL/hr
30 mL	309650	n/a	SS-30L (35 mL total)	n/a	0.1 to 650 mL/hr
35 mL	n/a	SoftPack: 1183500777 Rigid Pack: 8881535762	SS-30L (35 mL total)	n/a	0.1 to 650 mL/hr
50 mL	n/a	n/a	n/a	30602N	0.1 to 999 mL/hr
60 mL	309653	SoftPack: 1186000777 Rigid Pack: 8881560125	SS-60L	n/a	0.1 to 999 mL/hr

\*Syringes manufactured only in the United States

For product support, please contact Customer Advocacy at **888.812.3266** or email [customerfeedback@carefusion.com](mailto:customerfeedback@carefusion.com).

For technical support, please contact Instrument Technical Support at **888.812.3648**.

For product orders, please contact Customer Order Management at **800.482.4822**.

© 2012 CareFusion Corporation or one of its subsidiaries. All rights reserved. Alaris and IVAC are trademarks or registered trademarks of CareFusion Corporation or one of its subsidiaries. All other trademarks are property of their respective owners. IF1270-01 (1211/300)

CareFusion  
San Diego, CA

[carefusion.com](http://carefusion.com)

